
C|S

ELEMENTS

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ELEMENTS ANKLE SUPPORT

Carefully read this instruction sheet before use!

Indications / Intended Use

This product is intended for treatment of mild and temporary conditions such as pain and ankle sprains. The support can be used with or without the elastic straps. It gives a light compression and keeps the ankle joint warm. The movement of the ankle joint is limited from the elastic straps. This support is only intended for one person use.

Contraindications

Severe edema, Severe swelling, Neuropathy, Foot problems arising from diabetes, leg ulcers.

Adjustments

This product should only be altered by a medical professional person.

Warnings

Check that you have the correct size. Some conditions can result in swelling that can mean variations in size. If this is the first time you are using a support it is important to start slowly. Use the support for a short period of time the first day, then increase the use gradually. The support should not be used overnight. Do not use directly on open wounds. Use only on dry skin. Do not use lotion on your skin. Make sure the straps are not too tight. If any sign of skin irritation occurs, immediately remove the support. Check your skin on a daily basis.

Material Specifications

This product is made from: Polyamid, polyurethane and spandex. This product does not contain Latex, and is not flame retardant.




Washing Instruction

Close all Velcro prior to laundering. Machine wash in 40°C (104°F). Use a washing bag. Dry flat. Do not tumble dry or dry clean.

Who to Contact

If any problems occur, remove the support and contact your fitter or store where purchased.

Instructions for Use

1. Check that the support is complete and without damage.

2. Pull the sock on to the foot. Make sure it is comfortable.

3. The strap can be placed differently under the foot depending on how strong the limitation should be. If the straps are placed close to the heel the limitation will be less, and if it is placed further towards the toes it will be stronger. Be attentive that a too much limitation can be uncomfortable.

4. Pull the two straps around the ankle in a figure 8-formation. Secure in the front.

5. Fold the upper part of the sock over the straps for a sportier look.

6. Adjust for comfort.

7. If you need a stronger pull on only one side, both straps can be used on the outside (or inside).

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© Juli 2011