allard
Support for better life!

Instructions for Use

**Elements
Stretch Up**
Art. No. 34108

Instructions For Use

EN	ELEMENTS STRETCH UP	p. 4
SV	ELEMENTS STRETCH UP	p. 5
DA	ELEMENTS STRETCH UP	p. 6
FI	ELEMENTS STRETCH UP	p. 7
NO	ELEMENTS STRETCH UP	p. 8
DE	ELEMENTS STRETCH UP	p. 9
NL	ELEMENTS STRETCH UP	p. 10
IT	ELEMENTS STRETCH UP	p. 11
ES	ELEMENTS STRETCH UP	p. 12

Carefully read this instruction sheet before use!

Indications / Intended Use

This product is intended for treatment of Plantar Fasciitis. The sock should be used during night or at rest at any other time. This support is only intended for one person use. Single Patient-multiple use.

Contraindications

Circulation problems, open wounds on lower leg or foot, diabetes, skin problems or during pregnancy. Severe edema, Severe Swelling. Should not be used directly after any kind of surgery in leg or foot.

Adjustments

This product should only be altered by a medical professional person.

Warnings

Check that you have the correct size. Some conditions can result in swelling that can mean variations in size. If this is the first time you are using a support it is important to start slowly. Use the support for a short period of time the first day, then increase the use gradually. Do not use directly on open wounds. Use only on dry skin. Do not use lotion on your skin. Make sure the straps are

not too tight. If any sign of skin irritation occurs, immediately remove the support. Check your skin on a daily basis.

Material Specifications

This product is made from: 100% Spandex. This product does not contain Latex, and is not flame retardant.

Washing Instruction

Close all Velcro prior to laundering. Machine wash in 40°C (104°F). Use a washing bag. Dry flat. Do not tumble dry or dry clean. Wash with similar colours.

Disposal





Disposal in accordance with local legislation.

Information

Any serious incident that may occur in relation to the device should be reported to the manufacturer and the competent authority of the Member State in which the user and / or patient is established.

Who to Contact

If any problems occur, remove the support and contact your fitter or store where purchased.

Instructions for Use

1. Check that the support is complete and without damage.
2. Open the Velcro straps.
3. Pull the sock on to the leg. Make sure the surface is smooth without wrinkles.
4. Attach the Velcro around the leg under the knee.
5. Put the long strap through the loop and pull to stretch the sole of the foot without hurting. Close the Velcro. Start the treatment by pulling until the toes are approx. 1 cm from the floor. Increase the pull gradually - but not more the up til 2,5 cm difference between toes and floor.
6. Adjust for comfort.

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