## EN - INSTRUCTIONS FOR USE

### General Description

Manosplint® is a low temperature thermoplastic material used to fabricate splints, orthoses and rehabilitation aids.

### Precautions

To be used under the guidance of a qualified medical professional.

Any serious incident that has occurred in relation to the device should be reported to the manufacturer and the competent authority of the Member State in which the user and/or patient is established.

### Heating Methods, Temperature & Working Time

**Hot Water Method:** 140° - 162°F (60° - 72°C) for 1-2 minutes with a forming working time of 2 - 4 minutes depending on material thickness and the degree in which the material was stretched. Manosplint® Ohio and Manosplint® Ohio S turn translucent when heated. Manosplint® Ohio F should be placed in a sealed plastic bag to prevent water absorption into the fabric material.

**Dry Heat Oven Method:** 176° - 194°F (80° - 90°C) for 3-5 minutes with a forming working time of 2 - 4 minutes depending on material thickness and the degree in which the material was stretched. Manosplint® Ohio and Manosplint® Ohio S turn translucent when heated.

### Patterns & Cutting

A pattern can be drawn directly on the sheet material with a scratch awl. Thicker material may require scoring first with a utility knife and then bending at the score. Cut the pattern with a heavy duty shears or scissors. Finish trimming with a scissors.

### Symbols used

	Follow the instructions for use		CE marking of conformity		Medical device		Manufacturer		Manufacture batch number		Device reference		Temperature limits during storage and transport
---	---------------------------------	---	--------------------------	---	----------------	---	--------------	---	--------------------------	---	------------------	---	---

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Note, warming the material to soften it slightly (not heating all the way) will aid in the ease of cutting and form a soft sealed edge.

### Forming the Splint

Utilize desired heating method. Note Moderate resistive materials require smooth strokes while maximum resistive material will require a firmer hand. Manosplint® Ohio S and Manosplint® Carolina are uncoated material "sticky" this will aid in application by keeping the material in place. It will also allow for a permanent bond. Coated materials Manosplint® Arizona Manosplint® Ohio and Manosplint® Wisconsin can form a temporary bond by pinching the material together or a permanent bond if the coating is removed by scraping or with a solvent.

### Finishing & Adjusting

Manosplint® thermoplastic splints can be adjusted with a heat gun or dipping in hot water. Edges will remain smooth and closed. Splints maintain their shape while using low heat for fitting and adjustments. Manosplint® Ohio and Manosplint® Ohio S have 100% memory and if reheated can return to their original shape.

### Care

Clean with a mild soap and lukewarm water. Do not immerse fabricated splints in water greater than 123°F (51°C) as the splint will lose its shape.

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