Tro-Camp 21021

SVENSKA

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ENGLISH

Indications / Intended Use

Tro-Camp is a sacroiliac belt suitable for most figure types. It is intended to be used for hypermobile joints, and is suitable during pregnancy.

Contraindications

No found contraindications.

Adjustment

This product should only be altered by a professional orthotist.

Warnings

Do not pull the straps too tight, as this can cause bad circulation or skin irritation. If you feel any irritations or other problems, take it off immediately. Check skin on a daily basis.

Material specifications

Cotton textile, lined with cotton terrycloth.


Washing instruction

Machine washable in 40°C (104°F). Washing bag is recommended. Stretch to correct size when still wet.

Who to contact


If any problems occur please contact your orthotist or fitter.

Instructions for use


1. Undo the hooks and eyes on the front. Loosen up the buckles in the back and lengthen the straps. This makes the corset larger and easier to put on. The straps should still be positioned in the buckles.

2. Make sure that terrycloth is towards the skin.


3. It is helpful to be standing in front of a mirror when putting on the corset. Attach the hook and eyes in the front. Make sure the corset is positioned across the hip joint.


4. The diagonal strip of hooks and eyes is intended to accommodate the abdomen as it expands during later months of pregnancy.

5. Pull the straps in the back so that a comfortable support is achieved. NOTE! Best effect is achieved in a laying down position.



6. To remove the corset, undo the buckles in the back by turning them backwards and release the tension and the lacing can be opened.

7. If the corset tends to migrate, attach the groin straps in the push buttons in the front and the buckles on the back.

Item No.	Hip circumference	Height front	Height back
21021 0075	75 cm	10 cm	8 cm
21021 0080	80 cm	10 cm	8 cm
21021 0085	85 cm	10 cm	8 cm
21021 0090	90 cm	10 cm	8 cm
21021 0095	95 cm	10 cm	8 cm
21021 0100	100 cm	10 cm	8 cm
21021 0105	105 cm	10 cm	8 cm
21021 0110	110 cm	10 cm	8 cm

SVERIGE	DENMARK	SUOMI	NORGE	UK/IRELAND	USA	INTERNATIONAL
Tel 042-25 27 01	43 96 66 99	09 350 76 30	23 23 31 20	+353 42 932 8177	888-678 6548	+46 42 25 27 01
Fax 042-25 27 25	43 43 22 66	09 350 76 388	23 23 31 21	+353 42 932 8182	800-289 0809	+46 42 25 27 25
info@camp.se	info@camp.dk	info@camp.fi	info@camp.no	info@allarduk.co.uk	info@allardusa.com	info@allardint.com

Support For Better Life • www.allardsupport.com