









Reinforcements 39320 Trousers


elements
body



Three reinforcements are included in the price. If the box for Standard Panel is ticked this will count as one, and another two can be added.

Reinforcements 39320 Trousers


elements
body



IRP-R

Internal Rotation Panel


To assist internal rotation at
RIGHT hip



IRP-L

Internal Rotation Panel


To assist internal rotation at
LEFT hip



ERP-R

External Rotation Panel

To assist external rotation at
RIGHT hip




ERP-L

External Rotation Panel

To assist external rotation at
LEFT hip

Reinforcements 39320 Trousers

elements
body




IRPL-R

Long Internal Rotation Panel

To assist internal rotation at
RIGHT hip

Only for Long Legs




IRPL-L

Long Internal Rotation Panel

To assist internal rotation at
LEFT hip

Only for Long Legs




ERPL-R

Long External Rotation Panel

To assist external rotation at
RIGHT hip

Only for Long Legs



ERPL-L

Long External Rotation Panel



To assist external rotation at
LEFT hip

Only for Long Legs

Reinforcements 39320 Trousers



elements
body





PHE-L
Posterior Hip Extension Left
To assist hip extension
on LEFT side

PHE-R
Posterior Hip Extension
Right
To assist hip extension

AKE-L
Anterior Knee Extension
Left
To assist knee extension
on LEFT side

AKE-R
Anterior Knee Extension
Right
To assist knee extension

PKF-L
Posterior Knee Flexion Left
To assist knee flexion
on LEFT side

PKF-R
Posterior Knee Flexion Right
To assist knee flexion
on RIGHT side
Long Legs and ¾ Legs only

Three reinforcements are included in the price. If the box for Standard Panel is ticked this will count as one, and another two can be added.