


**39300, 39301 Vest
39305, 39306 Body
Reinforcements**

elements
body



SP
Standard Panel

To assist hip external rotation and back extension

LTPR
Lateral Trunk Panel

To assist RIGHT lateral flexion

LP-X
Lumbar Panel Extended

To assist postural control

SPL
Scoliosis Panels Left

To assist correction of right thoracic and LEFT lumbar curves

CENTER
Right thoracic
Left Lumbar

Please ensure that the correct reinforcement is selected.

LP
Lumbar Panel

To assist postural control

SPR
Scoliosis Panels Right

To assist correction of left thoracic and RIGHT lumbar curves

CENTER
Left thoracic
Right Lumbar


Please ensure that the correct reinforcement is selected.

LTPL
Lateral Trunk Panel

To assist LEFT lateral flexion

**39300, 39301 Vest
39305, 39306 Body
Reinforcements**

elements
body



SPL-C3
Scoliosis Panels Left
To assist correction of *LEFT* C-curve and right trunk rotation

SPR-C3
Scoliosis Panels Right
To assist correction of *RIGHT* C-curve and left trunk rotation

SPL-C2
Scoliosis Panels Left
To assist correction of *LEFT* C-curve without right trunk rotation.



SPR-C2
Scoliosis Panels Right
To assist correction of *RIGHT* C-curve without left trunk rotation.

SPL-C1
Scoliosis Panels Left
To assist correction of *LEFT* C-curve with mild right trunk rotation

SPR-C1
Scoliosis Panels Right
To assist correction of *RIGHT* C-curve with mild left trunk rotation

39300, 39301 Vest 39305, 39306 Body Reinforcements

elements
body



SBEL


Shoulder Back Extension Left
To assist retraction of the left scapulae and back extension on the RIGHT side

SBER

Shoulder Back Extension Right
To assist retraction of the right scapulae and back extension on the LEFT side

SBEL

SBER




PBP

Posterior Back Panel

To assist trunk extension

Model:
Paediatric,
Male

Model:
Female



AAP

Anterior Abdominal Panel



To assist trunk flexion

For Body Only

Three reinforcements are included in the price. If the box for Standard Panel is ticked this will count as one. Another two can be added. If Standard Panel is not ticked, three reinforcements can be chosen.

3930I Vest 39306 Body Reinforcements


elements
body



AEP-L
Anterior Elbow Panel Left

*To resist elbow hyperextension
on left side*


Long Sleeves Only



AEP-R
Anterior Elbow Panel Right

*To resist elbow hyperextension
on right side*


Long Sleeves Only



PEP-L
Posterior Elbow Panel Left

*To resist elbow flexion on left
side*

Long Sleeves Only



PEP-R
Posterior Elbow Panel Right

*To resist elbow flexion on right
side*


Long Sleeves Only

Three reinforcements are included in the price. If the box for Standard Panel is ticked this will count as one. Another two can be added. If Standard Panel is not ticked, three reinforcements can be chosen.

CAMP[®]
SCANDINAVIA

3930I Vest 39306 Body Reinforcements

elements
body



AS-L
Anterior Supination Left

*To resist supination
on left side*


Long Sleeves Only



AS-R
Anterior Supination Right

*To resist supination
on right side*


Long Sleeves Only



PP-L
Posterior Pronation Left

*To resist pronation
on left side*

Long Sleeves Only



PP-R
Posterior Pronation Right

*To resist pronation
on right side*

Long Sleeves Only

Three reinforcements are included in the price. If the box for Standard Panel is ticked this will count as one. Another two can be added. If Standard Panel is not ticked, three reinforcements can be chosen.

CAMP[®]
SCANDINAVIA