


Spare parts and accessories

MILINE is the next generation of modular, high-performance, orthotic ankle components from Becker Orthopedic. Offering superior durability, function and application in thermoplastic and composite orthoses. High strength fatigue and corrosion resistant steel. Heavy duty options available and extendable service life.

MILINE™ Spare parts



MILINE™ Standard Spring Pack


CAMP Item.no	Description	Size	Pack Contains:
625570013	MILINE Standard Spring Pack	13 mm	(1) 140 Rate spring (1) 140 Motion limiter pin (1) Ball bearing
625660016	MILINE Standard Spring Pack	16 mm	(1) 600 Rate spring, (1) 600 Motion limiter pin (1) Ball bearing (1) 320 Rate spring (1) 320 Motion limiter pin (1) Ball bearing
625750020	MILINE Standard Spring Pack	20 mm	(1) 320 Rate spring (1) 320 Motion limiter pin (1) Ball bearing



MILINE™ Pivot Bushing

Pivot bushing for **all** MILINE ankle joints.


CAMP Item.no	Description	Size	
625510013	MILINE Pivot Bushing	13 mm	
625600016	MILINE Pivot Bushing	16 mm	
625690020	MILINE Pivot Bushing	20 mm	



MILINE™ Pivot Screw

Pivot screw for **all** MILINE ankle joints.

CAMP Item.no	Description	Size	
625520013	MILINE Pivot Screw	13 mm	
625610016	MILINE Pivot Screw	16 mm	
625700020	MILINE Pivot Screw	20 mm	



MILINE™ Adjustment Pack

Adjustment set screw for MILINE Double Action and Dorsiflexion Assist ankle joints.

CAMP Item.no	Description	Size	
625530013	MILINE 13mm Adjustment set screw pack	13 mm	
625620016	MILINE 16mm Adjustment set screw pack	16 mm	
625710020	MILINE 20 mm Adjustment set screw pack	20 mm	



MILINE™ Adjustment Pack

Internal spring and parts for MILINE Double Action and Dorsiflexion Assist ankle joints with Booster Springs.


CAMP Item.no	Description	Size	Pack Contains:
625540013	MILINE Booster Spring Pack	13 mm	(1) 320 rate spring (1) 600 rate spring (1) Spring spacer (1) Booster hat and pin assembly
625630016	MILINE Booster Spring Pack	16 mm	(1) 1800 rate spring (1) 320 rate spring (1) Booster hat and pin assembly (1) Booster 1800 replacement spacer
625720020	MILINE Booster Spring Pack	20 mm	(1) 1800 rate spring (1) 320 rate spring (1) Booster hat and pin assembly (1) Booster 1800 replacement spacer



MILINE™ Coupler Pin Pack

Coupler pin for MILINE Double Action ankle joints.


CAMP Item.no	Description	Size	
625550013	MILINE Coupler Pin Pack	13 mm	
625640016	MILINE Coupler Pin Pack	16 mm	
625730020	MILINE Coupler Pin Pack	20 mm	



MILINE™ Hardware Pack

Hardware pack for all MILINE ankle joints.

CAMP Item.no	Description	Size	Pack Contains:
625560013	MILINE Hardware Pack	13 mm	(1) Pivot bushing (1) Pivot screw (2) Attachment screws (1) Teflon grease (1) Thread Lock
625650016	MILINE Hardware Pack	16 mm	(1) Pivot bushing (1) Pivot screw (2) Attachment screws (1) Teflon grease (1) Thread Lock
625740020	MILINE Hardware Pack	20 mm	(1) Pivot bushing (1) Pivot screw (2) Attachment screws (1) Teflon grease (1) Thread Lock



MILINE™ Thrust Plate Pack

Thrust plates for adjusting M/L play for all MILINE ankle joints.

CAMP Item.no	Description	Size	Pack Contains
625590013	MILINE Thrust Plate Pack	13 mm	(2) .010" thrust plates (2) .015" thrust plates (2) .020" thrust plates (2) .025" thrust plates Thrust plate installation tool
625680016	MILINE Thrust Plate Pack	16 mm	(2) .010" thrust plates (2) .015" thrust plates (2) .020" thrust plates (2) .025" thrust plates (2) .030" thrust plates Thrust plate installation tool
625770020	MILINE Thrust Plate Pack	20 mm	(2) .010" thrust plates (2) .015" thrust plates (2) .020" thrust plates (2) .025" thrust plates Thrust plate installation tool




MILINE™ Motion Stop Pin Pack

Motion stop pin for all MILINE Double Action and Dorsiflexion Assist ankle joints.

CAMP Item.no	Description	Size	
625580013	MILINE Motion Stop Pin Pack	13 mm	
625670016	MILINE Motion Stop Pin Pack	16 mm	
625670016	MILINE Motion Stop Pin Pack	20 mm	


MILINE™ Accessories



MILINE™ Short Upright

MILINE Short Uprights are designed specifically for use in bilateral orthotic designs.

CAMP Item.no	Description	Size	
62529 0013	MILINE lamination upright	13 mm	
62530 0016	MILINE lamination upright	16 mm	
62531 0020	MILINE lamination upright	20 mm	



MILINE™ Rivet Stirrup

Rivet stirrup for all MILINE ankle joints. Heavy-duty 16mm option available for patients with profound coronal plane postural abnormalities, or for those with excessive internal, or external foot progression angle.

CAMP Item.no	Description	Size	
625320013	MILINE rivet stirrup	13 mm	
625330016	MILINE rivet stirrup	16 mm	
625340020	MILINE rivet stirrup	20 mm	
625350016	MILINE heavy duty rivet stirrup	16 mm	



MILINE™ Standard Action Rivet Stirrup

Rivet stirrups for MILINE Standard Action ankle joints are laser etched in 5° increments for optimal adjustment of motion stop angle. Heavy-duty 16mm option available for patients with profound coronal plane postural abnormalities, or for those with excessive internal, or external foot progression angle.

CAMP Item.no	Description	Size	
625360013	MILINE stand act rivet stirrup	13 mm	
625370016	MILINE stand act rivet stirrup	16 mm	
625380020	MILINE stand act rivet stirrup	20 mm	
625390016	MILINE heavy duty stand act rivet stirrup	16 mm	



MILINE™ Standard Action Y-Stirrup

Y-stirrups for MILINE Standard Action ankle joints are laser etched in 5° increments for optimal adjustment of motion stop angle. Heavy-duty 16mm option available for patients with profound coronal plane postural abnormalities, or for those with excessive internal, or external foot progression angle

CAMP Item.no	Description	Size	
625400013	MILINE stand action y-stirrup	13 mm	
625410016	MILINE stand action y-stirrup	16 mm	
625420020	MILINE stand action y-stirrup	20 mm	
625430016	MILINE heavy duty stand action y-stirrup	16 mm	



MILINE™ Y-Stirrup

Y-stirrups for MILINE Standard Action ankle joints are laser etched in 5° increments for optimal adjustment of motion stop angle. Heavy-duty 16mm option available for patients with profound coronal plane postural abnormalities, or for those with excessive internal, or external foot progression angle

CAMP Item.no	Description	Size	
625440013	MILINE y-stirrup	13 mm	
625450016	MILINE y-stirrup	16 mm	
625460020	MILINE y-stirrup	20 mm	
625470016	MILINE heavy duty y-stirrup	16 mm	



MILINE™ Unilateral Short Upright

MILINE Unilateral Short Uprights are designed specifically for use in unilateral composite orthotic designs.

CAMP Item.no	Description	Size	
625480013	MILINE unilateral short upright	13 mm	
625490016	MILINE unilateral short upright	16 mm	
625500020	MILINE unilateral short upright	20 mm	