

# Massage Form Shaper Lightweight

Skalprotes/delprotes gjord i mjukt lättviktssilikon med masserande baksida. Kanalerna på protesens insida ger en skonsam massage mot bröstkorgen och tillåter luft att cirkulera, vilket håller huden sval och låter den andas.

Modellen är något extra fyllig i mitten och tunnas gradvis ut mot mjuka, tunna kanter för ett naturligt utseende i bh:n. Den runda formen gör det möjligt att vinkla protesen åt båda hållen för bästa möjliga passform.

## Användningsområde

Fördelaktig för kvinnor med lymfödem. Protesen lämpar sig väl för aktivitet och sport.



Art.nr.	Modell	Storlek
113850001	Massage Form Shaper Lightweight	1
113850002	Massage Form Shaper Lightweight	2
113850003	Massage Form Shaper Lightweight	3
113850001	Massage Form Shaper Lightweight	4
113850005	Massage Form Shaper Lightweight	5
113850006	Massage Form Shaper Lightweight	6
113850007	Massage Form Shaper Lightweight	7
113850008	Massage Form Shaper Lightweight	8
113850009	Massage Form Shaper Lightweight	9
113850010	Massage Form Shaper Lightweight	10
113850011	Massage Form Shaper Lightweight	11
113850012	Massage Form Shaper Lightweight	12