

Tawny Massageform Super Soft

En symmetrisk bröstprotes i mörkare nyans som ger mjuk massage vid användning. Kanalerna på protesens baksida masserar bröstkorgens vägg och tillåter luft att cirkulera så att huden känns sval och kan andas. Protesen är bland annat fördelaktig för kvinnor med lymfödem. Den speciella baksidan gör att protesen sitter säkert mot bröstkorgen. Framsidan är gjord av en extra mjuk silikon som ger en naturlig profil med mjukt fall. Tawny är 20 % lättare än traditionella silikonproteser. Bäst effekt fås i en stödjande BH. Protesen passar bäst den mindre och medelfylliga bysten dvs. A och B-kupa.

Användningsområde

Bla. för kvinnor som ligger i riskzonen att få lymfödem eller har lymfödem. Protesen lämpar sig väl för aktivitet och sport.



ABC Breast Care

Art.nr.	Modell	Storlek
102760001	Tawny Massage Super Soft	1
102760002	Tawny Massage Super Soft	2
102760003	Tawny Massage Super Soft	3
102760004	Tawny Massage Super Soft	4
102760005	Tawny Massage Super Soft	5
102760006	Tawny Massage Super Soft	6
102760007	Tawny Massage Super Soft	7
102760008	Tawny Massage Super Soft	8
102760009	Tawny Massage Super Soft	9
102760010	Tawny Massage Super Soft	10
102760011	Tawny Massage Super Soft	11
102760012	Tawny Massage Super Soft	12
102760013	Tawny Massage Super Soft	13
102760014	Tawny Massage Super Soft	14